



## Hitachi TM3000 - SEM Exploration Activity


Use the SEM Operation Instructions located in Dropbox or in the folder in front of the SEM.



### SafeSense

SEM Exploration – Specimen Analysis Sheet

Christopher Barr  
Cory Bogas  
Victoria Christensen



Sample ID	Line 1: Specimen Description Line 2: Purpose for Analysis	Page number in Journal
A	PMMA	
	To test shock absorbing characteristics	
B	300k MW PEO	
	To test shock absorbing characteristics	
C	900k MW PEO	
	To test shock absorbing characteristics	
D	Waterlock	
	To test shock absorbing characteristics	
<b>NOTES:</b> <ul style="list-style-type: none"> <li>Flow rate when we spun 300K PEO was too high</li> <li>Waterlock absorbs water</li> <li>PMMA is used to make plexi-glass</li> <li>PEO is a biocompatible material</li> </ul>		<b>APPROVAL:</b> <hr/>

### OBSERVATIONS:

#### Specimen #1:

- Less dense than other fibers
- fuzzy

#### Specimen #2:

- Thin mat of fibers
- Water droplets visible on fibers (high flow rate)

#### Specimen #3:

- When placing the fibers on the stub, the tweezers pinched the nanofibers, creating a rip in them

#### Specimen #4:

- Opaque fibers with water absorbency properties



Specimen 1 A (PMMA), Point 1

**100X**

SafeSense 0000 2014/02/28 08:39 A.L D5.3 x100 1 mm  
A (PMMA)

**OBSERVATIONS:**

- Random orientation
- Large Fibers
- Thin fibers (thickness)

**500X**

SafeSense 0001 2014/02/28 08:42 A.L D5.2 x500 200 um  
A (PMMA)

**OBSERVATIONS:**

- Flat Fibers
- Have the potential to roll up
- Like blades of grass or dental floss

**1500X**

SafeSense 0002 2014/02/28 08:46 A.L D5.8 x1.5k 50 um  
A (PMMA)

**OBSERVATIONS:**

- So thin that they are transparent
- Melted by the electron beam
- Large gaps between fibers

Specimen 1 A (PMMA), Point 2

**100X**

SafeSense 0003 2014/02/28 08:51 A.L D5.9 x100 1 mm  
A (PMMA)

**OBSERVATIONS:**

- Broken fibers
- Possibly where the tweezers pinched them

**500X**

SafeSense 0004 2014/02/28 08:52 A.L D5.9 x500 200 um  
A (PMMA)

**OBSERVATIONS:**

- Fibers starting to melt
- Twisting and rolling up
- Look thicker than normal fibers (in diameter)
- Possibly strong

**1500X**

SafeSense 0005 2014/02/28 08:56 A.L D5.9 x1.5k 50 um  
A (PMMA)

**OBSERVATIONS:**

- Looks like nano-scale tubes or straws
- Since air is in the tubes, it could have good shock absorbing qualities

**CONCLUSION:**

It is easy to visualize that this polymer is used in plexi-glass because it looks strong. Definitely want to attempt to spin this in bulk and drop test on a thick mat to see how well it absorbs shock. The pockets of air lead us to believe that it could be useful for our research.



Specimen 2 B (300k PEO), Point 1

**100X**

SafeSense 0000 2014/02/28 11:56 A.L D6.4 x100 1 mm  
B (PEO 300K)

**OBSERVATIONS:**

- Water droplets are visible
- More dots than fibers
- Individual fibers not as visible

**500X**

SafeSense 0001 2014/02/28 12:01 A.L D6.4 x500 200 um  
B (PEO 300K)

**OBSERVATIONS:**

- Cloudy
- Dots
- Black lines could be fibers but is unclear

**1500X**

SafeSense 0002 2014/02/28 12:03 A.L D6.4 x1.5k 50 um  
B (PEO 300K)

**OBSERVATIONS:**

- Thin fibers connected to drops of PEO polymer
- No fibers because of low MW

Specimen 2 B (300k PEO), Point 2

**100X**

SafeSense 0003 2014/02/28 12:06 A.L D6.1 x100 1 mm  
B (PEO 300K)

**OBSERVATIONS:**

- Could be a surface scratch on the fibers
- Or a big fiber?

**500X**

SafeSense 0004 2014/02/28 12:09 A.L D6.1 x500 200 um  
B (PEO 300K)

**OBSERVATIONS:**

- Still see droplets of PEO
- Looks like stitches or an x-ray of a spine

**1500X**

SafeSense 0005 2014/02/28 12:11 A.L D6.1 x1.5k 50 um  
B (PEO 300K)

**OBSERVATIONS:**

- Droplets of water?
- Too blurry to make out any fibers

**CONCLUSION:**  
This polymer would not be helpful for our research because of the low density and virtually no pockets of air. Without any fibers, we would not be able to analyze how they deform after impact.





Specimen 1 C (900k PEO), Point 1

**100X**

SafeSense 0000 2014/03/03 09:43 A.L D5.4 x100 1 mm  
C (PEO 900K)

**OBSERVATIONS:**

- Looks like aligned fibers
- Cloudy in some spots

**500X**

SafeSense 0001 2014/03/03 09:45 A.L D5.7 x500 200 um  
C (PEO 900K)

**OBSERVATIONS:**

- Fibers are not aligned
- Visible fibers due to a larger MW

**1500X**

SafeSense 0002 2014/03/03 09:47 A.L D5.7 x1.5k 50 um  
C (PEO 900K)

**OBSERVATIONS:**

- Thin fibers
- Not as many air pockets as PMMA
- Dried droplets of PEO or water

Specimen 1 C (900k PEO), Point 2

**100X**

SafeSense 0003 2014/03/03 09:51 A.L D5.3 x100 1 mm  
C (PEO 900K)

**OBSERVATIONS:**

- Spot where tweezers pinched the fibers
- Wrinkled surface is visible
- Can see aluminum foil is visible underneath

**500X**

SafeSense 0004 2014/03/03 09:52 A.L D5.4 x500 200 um  
C (PEO 900K)

**OBSERVATIONS:**

- Fibers bunched together after they were broken
- Does not look like traditional fibers
- Looks like many microorganisms

**1500X**

SafeSense 0005 2014/03/03 09:54 A.L D5.4 x1.5k 50 um  
C (PEO 900K)

**OBSERVATIONS:**

- Tightly packed fibers
- Started to melt under electron beam
- Looks webbed


**CONCLUSION:**

This polymer is interesting when it bunches up and deforms, so it would be beneficial to test it under our drop test apparatus to see how it reacts to impact. It has less pockets of air, so we expect it to work better than the 300K PEO but not as well as the PMMA.



Specimen 2 D (Waterlock), Point 1

**100X**




SafeSense 0000 2014/03/03 10:02 A.L D5.6 x100 1 mm  
D (Waterlock)

**OBSERVATIONS:**

- Small but visible fibers
- Small pockets of air
- Randomly assembled

**500X**




SafeSense 0001 2014/03/03 10:05 A.L D5.6 x500 200 um  
D (Waterlock)

**OBSERVATIONS:**

- Dense (more fibers)
- Look similar to PMMA

**1500X**




SafeSense 0002 2014/03/03 10:07 A.L D5.6 x1.5k 50 um  
D (Waterlock)

**OBSERVATIONS:**

- Transparent fibers
- Wavy
- Seem to be easy to manipulate/bend/squish

Specimen 2 D (Waterlock), Point 2

**100X**




SafeSense 0003 2014/03/03 10:11 A.L D5.4 x100 1 mm  
D (Waterlock)

**OBSERVATIONS:**

- Could be where water fell on the fibers
- Bunch together in two spots

**500X**




SafeSense 0004 2014/03/03 10:13 A.L D5.4 x500 200 um  
D (Waterlock)

**OBSERVATIONS:**

- Fibers branching off
- Water could have been absorbed which led them to be bunched together

**1500X**



SafeSense 0005 2014/03/03 10:14 A.L D5.4 x1.5k 50 um  
D (Waterlock)

**OBSERVATIONS:**

- Fibers become thicker in diameter the closer they get to the droplet
- Surface of the droplet has folds
- Maybe fibers become one when they absorb water?

**CONCLUSION:**

Since there are more of these fibers than any other ones that we analyzed, we would expect that these would work the best with impact. However, if we were to put this material in a helmet, sweat would absorb in it and could possibly change the characteristics of the fibers.